This is a qualitative, non-automated test intended for use by medical professionals. HemoTypeSC is intended to screen individuals of any age, including newborns, for abnormal hemoglobin variants HbS and HbC, which may indicate predisposition to sickle cell disease.

INSTRUCTIONS FOR USE

IVD HemoTypeSC™ is a rapid test kit intended to determine the presence of hemoglobin A, S and C in whole blood (optionally EDTA blood) & capillary blood. The kit is intended for *in vitro* diagnostic use by health professionals. This test is a competitive lateral flow immunoassay incorporating monoclonal antibodies for determination of the presence of hemoglobins A, S, and C. HemoTypeSC™ provides rapid detection of hemoglobin phenotypes HbAA, HbSS, HbSC, HbCC, HbAS, and HbAC.

- HemoTypeSC™ results are not affected by fetal hemoglobin. They may be obtained at any clinically appropriate time following birth.
- Phenotype HbS/α-thalassemia will produce the same HemoTypeSC™ result as HbSS.
- Other hemoglobin variants (e.g., HbD and HbE) will produce the same HemoTypeSC™ result as HbA.

Observe applicable biohazard safety procedures when handling Test Strips and human blood. Do not reuse any test strips or blood sampling devices. Some immobilized Test Strip reagents are derived from human sources - refer to local and national guidelines for biohazard safety handling procedures.

The HemoTypeSC™ Test Kit Includes:

- Foil Pouch containing the following:
 - Vial of 50 single-use Test Strips
 - Vial of 50 single-use Blood Sampling Devices
– vial may contain extra
 - Three reusable Dropper Pipettes
- Instructions for Use

Materials Not Provided:

- Water – use drinking water, **NOT** saline solutions
- Timer
- * Lancing Device and Lancets
- * Test Vials (1.5 - 5 mL tubes or vials may be used)
- * Rack for holding Test Vials
* Available as part of the HemoTypeSC™ Accessory

PROCEDURE

- Using Dropper Pipette, add six (6) drops of water (approximately 250 microliters) to Test Vial. Place Test Vial in a compatible rack.
- Open Vial of Blood Sampling Devices, remove one Blood Sampling Device, and reclose Vial. Obtain blood sample – a small drop is sufficient. Touch the white pad of the Blood Sampling Device to blood sample, until the white pad absorbs the blood droplet. In case of using whole blood: Do not immerse the blood sampling device into the tube of blood. **Ensure that the entire white pad has turned red.**
- Insert Blood Sampling Device into Test Vial water and swirl to mix.
 - Sufficient swirling is essential** for blood to be properly transferred into Test Vial.
 - Check visually to ensure that water has become **pink or light-red** in color.
 - Leave Blood Sampling Device in the Test Vial after swirling.
- Open Vial of Test Strips, remove one Test Strip, and reclose Vial. Insert HemoTypeSC™ Test Strip into Test Vial with arrows pointing **DOWN**.
- Wait 10 minutes.
- Take HemoTypeSC™ Test Strip out of the Test Vial and read results. Compare Test Strip to Results Chart on reverse side of this document for reference.

HemoTypeSC™ can be adapted for use with dried blood spots (DBS) if stored at 10-25°C and away from moisture for up to 5 months or at -20-0°C and away from moisture for up to 2 years. **The DBS should be cut to be 3/16" (4.76 mm) diameter circle and must be completely red with blood.**

PROCEDURE (DRIED BLOOD SPOTS)

- Using Dropper Pipette, add six (6) drops of water (approximately 250 microliters) to Test Vial. Place Test Vial in a compatible rack.
- Cut out DBS to be 3/16" (4.76mm) in diameter.
- Completely submerge the DBS under the water inside the Test Vial.
- Wait 10 minutes.
- Open Vial of Test Strips, remove one Test Strip, and reclose Vial. Insert HemoTypeSC™ Test Strip into Test Vial with arrows pointing **DOWN**.
- Wait 10 minutes.
- Take HemoTypeSC™ Test Strip out of the Test Vial and read results. Compare Test Strip to Results Chart on reverse side of this document for reference.

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